

KS Implant

The KS System boasts a distinctive internal design that fortifies the strength of the implant neck's coronal aspect, a critical area prone to fracture. This unique design feature ensures increased resilience, offering enhanced protection against potential fractures and reinforcing the implant's overall longevity.



OPTIMIZED TREATMENT OUTCOMES

Facilitates precise implant placement and prosthetic restoration, contributing to improved aesthetic and functional outcomes for patients.



SIMPLIFIED INVENTORY

By utilizing a single platform, inventory management is streamlined, reducing complexity and costs associated with stocking multiple implant systems.




ENHANCED CLINICAL FLEXIBILITY

Provides clinicians with the flexibility to choose the appropriate implant size and prosthetic solution based on patient-specific needs, without compromising on compatibility or performance.




Product Details



Abutment can be fastened with only one hand thanks to the abutment holding system



Convenient single platform for surgery / prosthesis which reduces the inventory burden



Enhanced strength by 2.4 times from increased wall thickness, 15° morse taper, and deeper connctions