

SPECIFICALLY DESIGNED FOR PLACING TAPERED IMPLANTS

OSSTEM[®]
IMPLANT

Taper KIT

The Taper Kit is designed for **tapered implants**, enabling efficient **bone-level placement** in just a few steps. Clinicians benefit from **maximum convenience**, while patients benefit from reduced chair time.



EXCELLENT DRILLING POWER

Maximizing drilling power with distinguished drill design, ultimately avoiding intermediate drilling leading to smoother and faster drilling results.




SAFE AND STABLE IMPLANTATION

A unique drill design that skips intermediate drilling, resulting in a smooth yet quick drill process.




PREVENTS BONE HEATING

Excellent drilling power leads to skipping intermediate drilling which prevents bone heating by maintaining low drilling temperature.




Product Comparison (vs Straight Drill)



123 straight Drills

VS



Taper Drills

Taper KIT (OTSK)

Top panel components:

Torque Wrench TW30B



Depth Gauge OSDG

